

---

## **DAV VII Turbo Oral Aerosol Evacuation System**

### **Only \$995.00**

For Mercury Vapor, Ozone, BPA dust, abrasion powder, laser smoke, VOCs, airborne viruses...

**Disposable, intake liners for moisture control**

**Internal stiffener arm for extended non contact hose support**

**5 ft, 3" diameter hands free intake hose with 5" laminar flow opening**

**4 stage filtration for moisture, particulate, odor and vapor control**

**Compact, only 12" X 12" X 15", lightweight, only 35 pounds**

**Rugged powder coated scratch resistant metal housing**

**Powerful yet quiet 770 cfm, 68 db, 110v, 60hz vacuum**

**Casters for multiple operatory mobility**

**Pre-filter for particulate capture down to 10 micron**

**HEPA filter for particulate capture down to 0.3 micron**

**Granulated activated carbon filter for odors and vapors**

**Limited Lifetime Warranty**

