



**Contact Us Today**

**877-336-8247**

**[info@dentairvac.com](mailto:info@dentairvac.com)**

## **DAV VII Turbo Oral Aerosol Evacuation System**

**For Mercury Vapor, Ozone, BPA dust, abrasion powder, laser smoke, VOCs, airborne viruses...**

**Disposable, intake liners for moisture control**

**Internal stiffener arm for extended non contact hose support**

**5 ft, 3" diameter hands free intake hose with 5" laminar flow opening**

**4 stage filtration for moisture, particulate, odor and vapor control**

**Compact, only 12" X 12" X 15", lightweight, only 35 pounds**

**Rugged powder coated scratch resistant metal housing**

**Powerful yet quiet 770 cfm, 68 dbL, 110v, 60hz vacuum**

**Castered for multiple operatory mobility**

**Prefilter for particulate capture down to 10 micron**

**HEPA filter for particulate capture down to 0.3 micron**

**Granulated activated carbon filter for odors and vapors**

**Limited Lifetime Warranty**

